




Separation Analyser LUMiReader® X-Ray





Illuminating even the most concentrated emulsions & suspensions.

-  Pigment segregation
-  Real-time separation stability
-  Concentration gradients
-  Sediment & cake structure
-  (In-)Homogeneity of dispersions, powder mixtures & solid matrices
-  Phase separation
-  Direct instrumental method for shelf life
-  Packing density/Solidosity
-  Foam structure

LUMiReader® X-Ray using



1 X-RAY SOURCE 2 SAMPLE 3 SENSOR



- 1 X-RAY SOURCE
- 2 SAMPLE
- 3 SENSOR
- 4 5 6 7 KINETICS OF TRANSMISSION

The LUMiReader® X-Ray is the first dispersion analyser designed for one purpose only: to study dispersability, stability, separation, and consolidation phenomena for completely transparent to completely opaque emulsions or suspensions. For the first time ever, illuminate your sample instantaneously from top to bottom. Solve your most challenging dispersion problems with complete insight. Go places light cannot.

Our patented solution combines x-ray vision with the proven STEP-Technology®, permitting highest spatial resolution, short sampling intervals and powerful detection technique. The LUMiReader® X-Ray creates monochrome and parallel x-rays ¹ with the help of a special crystal, which transmits the entire sample cell ² of 20 mm height. More than 1600 detectors record the transmitted beam, giving an unprecedented resolution ³. Instantaneous transmission profiles across your whole sample are converted into stability and separation rates. Particle concentration, sediment packing density can also be calculated based on the extinction profiles.

The obtained Space- and Time-resolved Extinction Profiles allow you to measure and understand particle stability, phase separation, and sediment consolidation in systems that where up to now, unknown frontiers of science and technology ⁴⁻⁷.

The LUMiReader® X-Ray allows you to optimize your formulation, measure its stability, and perform accurate shelf life tests, regardless of the shape or concentration of your dispersed particles and droplets. In-situ, real-time, non-invasive and non-destructive.

Typical applications are cosmetics, pharmaceuticals, paint & pigments, construction materials and fillers, as well as mining, ceramics and petrol industry dealing with complex emulsions, slurries and sludges.

SEPView®



your window to dispersion analysis

- ▶ Windows-based with Ribbon User Interface
- ▶ Plug & play, pack & go
- ▶ Individual user customization
- ▶ Full SOP concept (Creation, capture, data analysis)
- ▶ 8 different tools to understand (quantify) even the most complicated dispersion:

Time lapse measurement replay

Dispersion fingerprint

Instability index

Clarification


Phase separation

Sedimentation and creaming velocities

First derivative of integral transmission

Linear X-ray attenuation coefficient

- ▶ Analysis templates
- ▶ Zoom in and out of any area of interest
- ▶ Comparison of old and new measurements
- ▶ Comprehensive database security and full audit log
- ▶ Complies with 21 CFR Part 11



Safety features


What we have done for your safety

A full protection security circle certified by PTB – the national metrology institute providing scientific and technical services and Federal Office for Radiation Protection (BfS) protects the user. So anyone can employ the instrument free of risk and without worries about the x-ray radiation. No specific training course about x-ray instrumentation, nor any special operation licenses are required by the user.


A thick lead layer shields the environment from the x-ray exposure inside the apparatus. Because of this and various other safety designs, there is no detectable energy dose outside the instrument. State of the art interlocks sentinelling the sample lid position and the x-ray tube shutter action. Opening the lid during a measurement causes an immediate shut down of the shutter and the x-ray source, triggered by the sentinel circuit independent from software and user.

Applications & Samples


Phase separation of MoS₂-Graphit suspension




Phase and sediment formation of polymer dispersion





Visual appearance and x-ray intensity profiles of oil sands after separation



Packing density gradients of lime after centrifugal consolidation



X-ray phase detection of ongoing clarification and sediment growth (left) within even opaque systems (right)



Benefits

extremely
high
concentrations

- ▶ In-situ analysis of transparent as well as opaque dispersions
- ▶ No dilution of emulsions or suspensions
- ▶ Signal does not depend on particle shape
- ▶ See/understand complex dispersion behaviour
- ▶ Study the various instability mechanisms
- ▶ High resolution of phase separation of multicomponent systems
- ▶ Detect concentration gradients within phases and sediment
- ▶ Determine mean and space resolved packing sediment densities
- ▶ Endless monitoring of sample behaviour for long-time storage information
- ▶ Use any continuous phase of dispersing agent
- ▶ Real time, non-invasive and non destructive
- ▶ High-end analyser for QC, process monitoring and R & D



Application Areas

Abrasives
Batteries
Carbon black
Catalysts
Coated particles
Cosmetics
Ceramics
Construction materials
Lubricants
Magnetic particles

Microemulsions
Metals
Mineral powders
Nanosuspensions
Oil sands
Paints
Pigments
Rigid foams
Silicone emulsions
Solid electrolytes



Specifications

Measurement principle	Monoenergetic X-ray attenuation
Phase separation	High concentrated dispersions (transparent or opaque)
In situ sediment analysis	Packing density & structure
Stability analysis	From seconds to days or weeks
Consolidation measurement	Also in combination with LUMiFuge & LUMiSizer
Conformity	ISO/TR 13097; ISO/TR 18811; ISO 18747-1, CFR 21 Part 11

Samples	Suspensions, Emulsions, Suspo-Emulsions, Sludges, Slurries, Foams & Powders
Channels	1 sample
Volume	0.3 ml to 1.6 ml
Concentration	Up to 100 Vol%
Particle	any shape, from nano to microscale, no density restriction

Source	Monoenergetic X-ray, 17.48 keV, max 20 W at 40 kV, air cooled
Monochromator	Graphite
Disturbance free	No moving parts
Dimensions (WxHxD), Weight	47 x 24 x 44 cm ³ , 25 kg
Power supply	24 V, 220 W, Adapter (100 V to 240 V) included
Safety	Fully radiation protected system; Radiation < 1 µSv/h (BfS 03/13 V RöV)
Radiation control requirements	None; instrument can be used anywhere



LUM GmbH, Berlin, Germany

Phone: +49 30 6780 60 30

E-Mail: info@lum-gmbh.de

Web: www.LUM-GmbH.com

www.lumireader-xr.com

www.dispersion-letters.com



The NEXT STEP in Dispersion Analysis & Materials Testing

Distributed by:

